

NPPTL COVID-19 Response: International Respirator Assessment

Manufacturer: Huizhou Huinuo Technology Co., Ltd.

Model Tested: 9501A

Date Tested: June 22, 2020

These findings pertain to the Huizhou Huinuo Technology Co., Ltd., model 9501A. The packaging for this product indicates that it meets GB2626-2006 (the Chinese standard for Respiratory Protective Equipment – Non-Powered Air-Purifying Particle Respirator).

Thirty respirators were submitted for evaluation. The respirators were sampled into groups of ten for evaluation. The samples were tested using a modified version of NIOSH Standard Test Procedure (STP) TEB-APR-STP-0059. This modified assessment plan can be found [here](#).

No certificate of approval was provided with the samples received; therefore, the authenticity of the claims cannot be validated.

The maximum and minimum filter efficiency was 98.10% and 96.23%, respectively. All thirty respirators measured more than 95%.

While the above-listed product classification has similar performance requirements to NIOSH-approved devices, NIOSH does not have knowledge about the sustained manufacturer quality system and product quality control for these products. NIOSH also does not have knowledge about the product's handling and exposures after leaving its manufacturer's control.

In addition, this product is an ear loop design. Currently, there are no NIOSH-approved products with ear loops; NIOSH-approved N95s have head bands. Furthermore, limited assessment of ear loop designs, indicate difficulty achieving a proper fit. While filter efficiency shows how well the filter media performs, users must ensure a proper fit is achieved.

This assessment is not a part of the NIOSH respirator approval process and will in no way lead to or preclude NIOSH approval through the official approval process. This assessment was developed as an assessment of the filter efficiency for those respirator's represented as certified by an international certification authority, other than NIOSH, to support the availability of respiratory protection to US healthcare workers due to the respirator shortage associated with COVID-19. Only particulate filter efficiency was assessed.

The results provided in this letter are specific to the subset of samples that were provided to NPPTL for evaluation.

These results will be used to update the CDC guidance for [Crisis Capacity Strategies \(during known shortages\)](#).

Evaluation of International Respirators

Test: Modified TEB-APR-STP-0059

Date Tested: June 22, 2020

Report Prepared: June 22, 2020

Manufacturer: Huizhou Huinuo Technology Co., Ltd.

Item Tested: 9501A (Sample Group 1 of 3)

Country of Certification: China (GB2626-2006)

Pictures have been added to the end of this report.

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
1	85	4.4	1.04	2.95	97.05
2	85	5.1	0.86	2.75	97.25
3	85	5.2	1.07	3.00	97.00
4	85	5.0	1.38	3.66	96.34
5	85	4.5	0.85	3.30	96.70
6	85	6.2	1.01	3.77	96.23
7	85	4.7	0.73	2.13	97.87
8	85	5.5	0.41	2.01	97.99
9	85	5.7	1.11	3.40	96.60
10	85	7.1	2.30	3.11	96.89
Minimum Filter Efficiency: 96.23			Maximum Filter Efficiency: 97.99		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

NPPTL COVID-19 Response: International Respirator Assessment

Test: Modified TEB-APR-STP-0059

Date Tested: June 22, 2020

Report Prepared: June 22, 2020

Manufacturer: Huizhou Huinuo Technology Co., Ltd.

Item Tested: 9501A (Sample Group 2 of 3)

Country of Certification: China (GB2626-2006)

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
11	85	5.4	0.92	2.65	97.35
12	85	6.0	0.65	2.56	97.44
13	85	5.0	0.83	2.81	97.19
14	85	4.4	0.88	2.69	97.31
15	85	5.7	0.84	2.69	97.31
16	85	4.8	0.64	2.06	97.94
17	85	5.7	0.64	1.90	98.10
18	85	6.0	0.54	3.00	97.00
19	85	4.9	0.69	2.05	97.95
20	85	4.6	0.90	2.14	97.86
Minimum Filter Efficiency: 97.00			Maximum Filter Efficiency: 98.10		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.

NPPTL COVID-19 Response: International Respirator Assessment

Test: Modified TEB-APR-STP-0059

Date Tested: June 22, 2020

Report Prepared: June 22, 2020

Manufacturer: Huizhou Huinuo Technology Co., Ltd.

Item Tested: 9501A (Sample Group 3 of 3)

Country of Certification: China (GB2626-2006)

Filter	Flow Rate (Lpm)	Initial Filter Resistance (mmH ₂ O)	Initial Percent Leakage (%)	Maximum Percent Leakage (%)	Filter Efficiency
21	85	5.3	0.82	3.38	96.62
22	85	4.6	0.98	2.52	97.48
23	85	6.1	0.74	2.29	97.71
24	85	6.0	1.11	3.10	96.90
25	85	4.9	1.10	2.40	97.60
26	85	5.5	1.26	3.65	96.35
27	85	6.2	0.88	3.39	96.61
28	85	4.9	0.97	2.82	97.18
29	85	4.6	0.84	2.72	97.28
30	85	5.3	0.68	2.72	97.28
Minimum Filter Efficiency: 96.35			Maximum Filter Efficiency: 97.71		

- The test method utilized in this assessment is not the NIOSH standard test procedure that is used for certification of respirators. Respirators assessed to this modified test plan do not meet the requirements of STP-0059, and therefore cannot be considered equivalent to N95 respirators that were tested to STP-0059.
- Respirators tested may not be representative of all respirators with the same certification mark. NIOSH has no control over suppliers and distributors of respirators certified by other national or international parties.
- This assessment is not a confirmation that it conforms with any or all of its specifications in accordance with its certification mark.
- This assessment was not a part of the NIOSH approval program. These results do not imply nor preclude a future approval through the NIOSH respirator approval program.



FACE MASK

GB2626-2006 KN95

9501A

Wearing Diagram of Mask



Hold the mask by hand, chin against mask, and cover the nose and mouth, ear strap Hanging on the ear root

Adjust the position of the mask, press the metal nose clip from the center of the Nose beam to both sides at the same time

Press it all the way to the bridge of the nose

Precautions for Use

The mask can be used for respiratory protection of some non oil particles (dust, pollen and other Kinds of productive dust and no volatile haze and other air pollutants)

This product can not be used in anoxic environment

This product can not be used for the protection of toxic gas and steam

If the mask is polluted, damaged, aged or invalid, it needs to be replaced

Children must use it under the guidance of adults. Please keep it in a safe place to avoid accidents!


Product Name: KN 95 Protective Mask
Model: 9501A
Expiry Date: 2-years (Original Packaging / Moisture Proof)

Manufacturer: Huizhou Huinuo Technology Co., Ltd
Address: No.26, Laohuling Group North Road, Xiaojinkou, Huicheng District, Huizhou City
www.hvn-china.com
TEL:86-0752-3195990



Made in China

NPPTL COVID-19 Response: International Respirator Assessment



NPPTL COVID-19 Response: International Respirator Assessment



NPPTL COVID-19 Response: International Respirator Assessment



NPPTL COVID-19 Response: International Respirator Assessment

